IsoMed 2000 product line

The IsoMed 2000 is a device for testing and training in static or dynamic mode.

These machines are successfully used in a wide range of biomechanical and therapeutic areas like:

Rehabilitation -Science and Research – or for Top athletes

This IsoMed 2000 system is equipped with a torque controlled servo motor system in its dynamometer head and a high precise telemetric torque-measurement unit.

This servo drive is currently the strongest on the market to guarantee a very stable movement performance with enough power reserves even when top athletes are going for broke on it.

Due to its construction and accessory adapters, belts or pads, it is possible to operate with all big joints and plans of the human body in a single isolated movement to evaluate the currently absolute or relative muscle stretch and power situation under static or dynamic conditions.

Muscle imbalances or mobilisation deficits often lead to increased joint degeneration and a much higher risk of injuries during physical activities in sports as well as by physical hard working peoples.



It is also possible to enhance training efficiency or modify therapy programs by analyzing the results of tests to retest comparison or continuous testing sequences in certain intervals.

In the fact of the dynamometer satellite concept and overhead monitoring combined with an unique automatic computer and motor controlled reposition system (*MEMOTRONIC*) for all seat or dynamometer adjustments it is not comparable to any other system designed to operate in highest accuracy and safety levels.



With the IsoMed2000 Basic equipment you can realise the following movements.

Ankle:

Knee


Hip

Hip - Flexion/Extension


Hip - Ab/Adduction


Shoulder

Shoulder - Flexion/Extension


Shoulder - Internal/External Rotation

Wrist

Wrist - Extension/Flexion

Wrist/Forearm - Pronation/Supination


Election of special joints


Rump Rotation

Power Lifting

Optional attachment: IsoMed2000 Back Module

IsoMed2000 Linear-module = “Legpress athletic” for working in a closed chain pattern

IsoMed2000 analyze program

Subjects

Quick Jump

Search Subject's Last Name for [] Search [] Last Results [] Create new Subject []

Subject overview

General Information View/Edit Details

Name: **REINHOLD FERSTL**
 Sex: **Male**
 Dominant Side: **right**

Subject details Cycles Evaluations

Create new Cycle []

Test Cycle View/Edit Details

General cycle data

Type: **Sport/Training cycle**
 Start: **10.11.2006 Ongoing**
 Cycle Rating: **Not available yet**
 Subject's weight: **80.0 kg**
 Subject's height: **170 cm**

Session blocks from 10.11.2006

Date/Time	Block Type	Joint	Test Mode	Plane	Block type
10.11.2006 21:03	Select/None	Knee left	isometric conc/conc	flexion	Select/None
10.11.2006 21:02	Select/None	Knee left	isometric conc/conc	flexion	Select/None
10.11.2006 20:59	Select/None	Knee right	isometric conc/conc	flexion	Select/None
10.11.2006 20:58	Select/None	Knee right	isometric conc/conc	flexion	Select/None

Evaluations

Gather data for the Evaluation

Below you can select the blocks you wish to include in the evaluation. You can adjust the selection by using the selection filters.

Evaluation parameters

Evaluation Type: **Left-Right Comparison**
 Sample group: **Intra one Subject**

Joint: **Knee**
 Test Mode: **isometric conc/conc**
 Plane: **flexion**
 Block type: **Select/None**

Clear Block Parameters []

Evaluation Container

Subject's name	External ID	Date of Birth	Age	
FERSTL, REINHOLD	-	-	-	Delete Subject []
10.11.2006 20:58	Knee right	60''s 60''s	2 Sets / 6 Reps	
10.11.2006 21:03	Knee left	60''s 60''s	2 Sets / 6 Reps	

Show/Hide selection filters []

Results of Subject filtering

Subject's name	External ID	Date of Birth	Age	
FERSTL, REINHOLD	-	-	-	
Test Cycle				
10.11.2006 21:02	Knee left	120''s 120''s	2 Sets / 6 Reps	Display [] Swap with FERSTL, REINHOLD: 10.11.2006 21:03 [] Swap []
10.11.2006 20:59	Knee right	120''s 120''s	2 Sets / 6 Reps	Display [] Swap with FERSTL, REINHOLD: 10.11.2006 20:58 [] Swap []

Subject Information

Enter the subject's information below. * = required field

General subject data

External ID []
 TIN []
 Last Name* **FERSTL**
 First Name* **REINHOLD**
 Date of Birth [] DD.MM.YYYY
 Sex **Male**
 Dominant Side **Right**
 Occupation []
 Insurance No []
 Statistical index: **No indexes to select from**
 Notes []

Contact data

Street []
 Code []
 Area []
 Phone []
 Fax []
 Email []

Custom Fields

Custom Field ONE []
 Custom Field TWO []
 Custom Field THREE []
 Custom Field FOUR []
 Custom Field FIVE []
 Custom Field SIX []

Save Subject Information [] Return to Subjects Overview []

analyze
Overview Subjects Evaluations Data Import Options

Quick Jump
Search Subjects: Last Name [] [] Search [] Last Results [] Create new Subject []

Subject overview

Isomed - Windows Internet Explorer

Short Overview

General Information
Date & time: 10.11.2006 21:03
Joint: Knee left
Test Mode: isokinetic conc/conc
Plan: flexion/extension
Motion Range: 10°-90°
Velocity: 60°/s, 60°/s
Break: 30 sec
Sets: 2
Repetitions: 6
Input Date: 14.11.2006

Editable Information
Block type: Select/None
Subject's weight: n.g. 60.5 kg
Subject's height: n.g. 178 cm
Reassign to: No matching cycles
Notes: []
Save Changes []

Additional Options
Show 3D-Overview []

Adjust data view
View Data & Charts for: whole block [] . show: isokinetic [] values in tabular view [] Change view []

Numeric Representation

	flexion	CoV	extension	CoV	flex/ext Ratio
Value	117.9 Nm		127 Nm		

analyze
Overview Subjects Evaluations Data Import Options

Quick Jump
Search Subjects: Last Name [] [] Search [] Last Results [] Create new Subject []

Subject overview

Isomed - Windows Internet Explorer

Fatigue in whole Block

Fatigue ratios in whole Block

Set	Flex (%)	Ext (%)	Mean Torque (Nm)
Set 1	40.1%	47.3%	84.0 Nm
Set 2	33.9%	37.4%	85.0 Nm

Mean ROM analysis flexion for whole Block
overall mean ROM 76.7 degrees

ROM Type	Value
acceleration ROM	5.9 deg
isokinetic ROM	60.7 deg
deceleration ROM	10.1 deg

Mean ROM analysis extension for whole Block
overall mean ROM 76.1 degrees

ROM Type	Value
acceleration ROM	5.9 deg
isokinetic ROM	60.9 deg
deceleration ROM	9.3 deg

analyze
Overview Subjects Evaluations Data Import Options

Quick Jump
Search Subjects: Last Name [] [] Search [] Last Results [] Create new Subject []

Subject overview

Isomed - Windows Internet Explorer

Parameter	6 Repts	0.11	6 Repts	0.11	30/32
Mean Power	88°	-	88°	-	50/50
Peak ROM	81°	-	81°	-	50/50

Peak Torque per Set in whole Block

Set	Flex (Nm)	Ext (Nm)
Set 1	112.2 Nm	137 Nm
Set 2	111.5 Nm	176 Nm

Peak Work per Set in whole Block

Set	Flex (Nm)	Ext (Nm)
Set 1	139.1	188.2
Set 2	118.1	187.1

Best Curves for Peak/Mean Torque in whole Block




flexion/extension ratio for whole Block

Parameter	Flex (%)	Ext (%)
Peak torque	38.3%	61.7%
Mean torque	37.4%	62.6%

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Berlin New paper BZ 2006

So is Gilberto agonizing

Sophisticated machine shall enable a faster comeback for the Brazilian

Gilberto works out hard for his comeback
Foto: ROGMANS

Swollen veins on his temples, sweat running down his brow – Gilberto works pretty hard at the Reha-Center “Da Vinci” (Katharinenstraße).. The Brazilian plugs away at his comeback – and one of the most modern and sophisticated trainings machines of the world support him on his task...

This machine named IsoMed 2000 is to build up the Brasilians’ muscles to Soccer League level again. “In Berlin we are the only Reha-Center owning such a machine. Normally it can be found in much bigger University hospitals. But due to this machine we can perform a very precise therapy.” expresses Thomas Sennewald, Physio in Chief.

Nearly every day Gilberto can be found flaying at this High-Tech machine. His dream as he speaks out: “At best I would be playing against Gladbach next Saturday again, but the doctor slows me down.”

And this is a rather good idea – just four weeks ago Gilberto broke his “Syndesmoseband”. Then the focus was on a seven week break. But the healing process went much better than first thought. It is next week that Gilberto shall train on the football ground again. He admires: “Gladbach might be much too early, that’s true. Even so the DFB-Cup against Stuttgart on Wednesday after this I won’t be on the field. But I am sincerely praying to be part of the Cottbus-Game on October the 28th.”

Until then Gilberto will agonize himself every day for the blitz-comeback...

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